**1. Intended use / Indication**

The **printo<sup>®</sup> GR-19.1 OA** resin is a light-curable polymerizable resin intended to be used in conjunction with extra-oral curing light equipment for the fabrication, by additive manufacturing, of dental objects such as mouthguards, nightguards, splints, repositioners and retainers.

**2. Contraindication**

**printo<sup>®</sup> GR-19.1 OA** is contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... for bruxism splints and splints with wall thickness ≤ 1.5 mm.
3. ... for denture bases.
4. ... for every application that is not part of the indication (see above).

**3. Patient Target Group**

Persons being treated in the context of a dental procedure.

**4. Intended Users**

Dentist, dental technicians

**5. Requirements**

**Software – information available from:**

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**Hardware (3D Printing) – information available from:**

ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

**Hardware (Post Curing) – information available from:**

pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

Additional specifications on [www.pro3dure.com](http://www.pro3dure.com)

**6. Material**

**printo<sup>®</sup> GR-19.1 OA** consists of functional (meth)acrylic resins, initiators, dyes and stabilisers.

**7. Geometric presetting**

Minimum wall thickness: 1.5 mm

**8. Material parameters**

Radiation penetration depth controllable by exposure time

50 µm

100 µm

**9. Manufacturing process (fig. 1-10)**

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts, until there are no residues of IPA or equivalent cleaning agent.
9. Post curing (4 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.

**10. Finishing processes**

Polishing

**11. Notice**

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. **Caution:** Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Avoid breathing dust/fume/gas/mist/vapours/spray. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed or concerned: Call a POISON CEN TER/doctor. If skin irritation occurs: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Collect spillage. Store locked up. Dispose of contents/container according to official regulations.


**12. Hazard statements**

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. Suspected of damaging the unborn child. Toxic to aquatic life with long lasting effects.

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exocad GmbH · Julius-Reiber-Straße 37, 64293 Darmstadt, Germany