



validierte 3D-Drucker | validated 3D printer

printo[®]
dent

		GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model	GR-13.3 model release	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary It	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex
ASIGA MAX		✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓
ASIGA MAX UV / MAX 2		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA PRO 2		✓	○	○	○	✓	✓	✓	✓	✓	○	○	✓	✓	✓	✓	○	○	○	○	○	○	○
ASIGA Ultra		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA PRO 4K		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
SprintRay Pro 95		○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
MiiCraft Ultra		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MiiCraft Prime		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MiiCraft Profession 4K		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MiiCraft Alpha		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Organical X1S		✓	○	○	○	✓	✓	✓	✓	✓	○	○	✓	✓	✓	✓	○	○	○	○	○	○	○
Organical X1N		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
phrozen Shuffle/Sonic		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓




validierte 3D-Drucker | validated 3D printer

printo[®]
dent

		GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model	GR-13.3 model release	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary It	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex
SISMA EVERES UNO		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Prusa MEDICAL ONE		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Anycubic Photon Mono X		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓
Microlay Versus		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓
Microlay Eve pro		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓
SHINING 3D AccuFab-L4D		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
SHINING 3D AccuFab-D1s		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓
Flashforge Focus 6K		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓
Structo ST-01		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Uniz SLASH PLUS		○	✓	✓	○	○	○	○	○	○	✓	✓	○	○	✓	○	✓	○	✓	○	○	○	○
UNIZ NBEE		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
NOBIL METAL Sinergia Print		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

✓ erhältlich / available ✓ in Bearbeitung / in progress ○ auf Anfrage / on request

pro3dure arbeitet kontinuierlich an der Entwicklung und Aktualisierung von Parametern für ausgewählte offene 3D-Drucksysteme./
pro3dure works continuously on the development and updating of parameters for selected open 3D printing systems.





validierte Nachhärtung | validated post curing

printo dent

	GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model	GR-14.2 denture release	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex	
pro3dure CD-1	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
pro3dure CD-2	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
NK-OPTIK Otoflash G171	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA Flash	○	○	✓	✓	✓	✓	✓	✓	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Kulzer cara Print LEDCure	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	○	○	○	○	○	○	○	○	○	○	○	○	○
SHINING 3D FabCure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	○	○	○	○	○	○	○	○
SprintRay ProCure	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
SprintRay ProCure 2	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
SprintRay NanoCure	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
SISMA EVERES CURO	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Prusa MEDICAL CW ONE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
WICKED ENGINEERING CUREbox	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

✓ erhältlich / available ✓ in Bearbeitung / in progress ○ auf Anfrage / on request



validierte Nachhärtung | validated post curing

printo dent

	GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model	GR-14.2 denture release	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex	
rapidshape RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Nobil Metal Dentafarm Photopol	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Meccatronicore BB Cure N	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓



validierte Reinigung | validated cleaning

printo dent

	GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model	GR-14.2 denture release	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex	
pro3dure CLD-1	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
SprintRay Pro Wash/Dry	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Kulzer cara Print Clean pro	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Prusa MEDICAL CW ONE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
rapidshape RS Wash	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Meccatronicore BB Wash	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

pro3dure arbeitet kontinuierlich an der Entwicklung und Aktualisierung von Parametern für ausgewählte Nachhärtungs- und Reinigungsgeräte./
pro3dure works continuously on the development and updating of parameters for selected post curing and cleaning devices.